



Saliva sampling made easy and valid!

Instructions on
how to use the app




Connecting to server



Distribute a QR code to the participant that contains the server address.

Only on first use of app.

Connecting to server



Ask participants to download the barcode scanner from Zxing Team.



Only on first use of app.


Connecting to server



The server address will appear and needs to be saved.

Only on first use of app.


Enter participant code



Distribute a code to the participants that can be used to associate the datasets.

Only on first use of app.

Enter participant code



Participants should type in their code and save it.


Only on first use of app.

Startscreen




Clicking here will provide participants with instructions on how to collect saliva samples and allow settings to be adjusted.

Menu




Instructions




The instructions are presented in the form of a slide show. Clicking Back will return participants to the startscreen.

Menu




Settings



Clicking here will
allow the adjustments
of settings.


Redirection to questionnaire



The app offers the option to combine saliva tracking with the assessment of other variables through an online questionnaire.

Only on first use of app.

Redirection to questionnaire



Distribute a QR code to the participant that contains the URL to the questionnaire if redirection is wanted.


Only on first use of app.

Startscreen




Clicking here will start the collection of the first saliva sample.

Collect sample




Once the cotton roll is in the mouth, participants can start the timer.

Collect sample




After 1min participants are asked to return the cotton roll and continue.

Photograph sample



Clicking here will
open the camera on
the phone.


Photograph sample



Ask participants to allow access to camera.

Only on first use of app.


Photograph sample



Ask participants to allow access file system.


Only on first use of app.

Photograph sample




Clicking here confirms the photograph. Label should be visible.

NO redirection to questionnaire



Clicking here will close the app.

Redirection to questionnaire



Clicking here will redirect participants to the questionnaire.